OsteoBiol[®]
by TecnoSS

Bone Grafting Materials

REGENERATION SCIENCE

INSPIRED BY NATURE

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Giuseppe Oliva MD
R&D Director
Tecnoss S.r.l.

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(1) Figueiredo M, Henriques J, Martins G, Guerra F, Judas F, Figueiredo H

Physicochemical characterization of biomaterials commonly used in dentistry as bone substitutes - comparison with human bone

Journal of Biomedical Materials Research Part B: Applied Biomaterials, 2010 Feb; 92(2):409-19

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
(2) Nannmark U, Sennerby L

The bone tissue responses to prehydrated and collagenated cortico-cancellous porcine bone grafts: a study in rabbit maxillary defects
Clinical Implant Dentistry and Related Research, 2008 dec;10(4):264-70

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Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy
Clinical Oral Investigations, 2018 Jan;22(1):505-513

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
(4) Rombouts C, Jeanneau C, Camilleri J, Laurent P, About I
Characterization and angiogenic potential of xenogeneic bone grafting materials: Role of periodontal ligament cells
Clinical Oral Investigations, 2018 Jan;22(1):505-513

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







(5) Brunelli G, Sollazzo V, Carinci F, Palmieri A, Girardi A, Monguzzi R
OsteoBiol® influences osteogenic differentiation of adipose derived stem cells
European Journal of Inflammation, 2011, Vol. 9, no. 3 (S), 103-107

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Wide diameter immediate post-extractive implants vs delayed placement of normal-diameter implants in preserved sockets in the molar region: 1-year post-loading outcome of a randomised controlled trial
European Journal of Oral Implantology, 2017;10(3):263-278

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
BIBLIOGRAFIA

- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC; 10(4):264-70
- (3) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (4) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42
- (7) MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVANI B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR
MATERIALS, 2017 SEP 8;10(9)
- (8) CHECCHI V, FELICE P, ZUCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278
- (9) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- (10) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (11) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRA-BONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (12) ROMBOUTS C, JEANNEAU C, CAMILLERI J, LAURENT P
ABOUT I CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENEIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

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BIBLIOGRAPHY

(1) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

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(3) GIULIANI A, IEZZI G, MAZZONI S, PIATELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVESTIG, 2017 2018 JAN;22(1):505-513

(4) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870

(5) IEZZI G, PIATELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44

(6) SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9

(7) BARONE A, BORGIA V, COVANI U, RICCI M, PIATELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

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
(9) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

(10) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68

(11) BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

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BIBLIOGRAFIA

- (1) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005 JUN;54(6):351-62
- (2) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006;2: 99-106
- (3) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (4) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (5) SANTAGATA M, GUARINIELLO L, TARTARO G
A MODIFIED EDENTULOUS EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N.:114-9
- (6) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

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mini-rialzo crestale
casi clinici a pag 82

BIBLIOGRAFIA

- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RESSIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67
- (5) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (6) LORENZON G, BUTTARELLO GM, CHESSA G
CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION
DENTISTRY, 2015, 5:2

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



RIALZO DI SENO PER VIA CRESTALE
mini rialzo con osteotomi
casi clinici a pag 82

BIBLIOGRAFIA

- (1) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- (2) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC; 77(12):1984-90
- (3) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL[®] INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107
- (4) CAKIR M, KARACA IR, AYSEGÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION
CLIN ORAL IMPLANTS RES, 2015 MAR-APR; 30(2):279-85
- (5) KOLMAS J, SZWAJA M, KOLODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENIC BIOMATERIALS USED AS BONE SUBSTITUTES
J PHARM BIOMED ANAL, 2012 MAR 5; 61:136-41
- (6) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U
CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL
J CLIN PERIODONTOL, 2017 FEB; 44(2):204-214
- (7) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006, 2: 99-106
- (8) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2015 DEC; 17(6):1114-26
- (9) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL
CLIN IMPLANT DENT RELAT RES, 2015 DEC; 17(6):1114-26
- (10) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC; 23(12):1369-76
- (11) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18

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BIBLIOGRAFIA

(1) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870

(2) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR-APR;33(2):159-66

(3) FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
IMPLANT DENT, 2012 JUN;21(3):190-5

(4) FELICE P, BARAUSSE C, BARONE A, ZUCHELLI G, PIATELLI M, PISTILLI R, IPPOLITO DR, SIMION M
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

(5) BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

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Packaging

20x15x5 mm, 20x10x5 mm

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
2506/14 | 20x15x5 mm | Soft | Suino curvo
2506/17 | 20x10x5 mm | Norm | Suino curvo

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






- **A)** Barone A et al. - **Flap versus flapless procedure for ridge preservation in alveolar extraction sockets: a histological evaluation in a randomized clinical trial**
Clinical Oral Implants Research, 2015 Jul;26(7):806-13
- **B)** Giuliani A et al. - **Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy**
Clinical Oral Investigations, 2018 Jan;22(1):505-513
- **C)** Barone A et al. - **Clinical and histological changes after ridge preservation with two xenografts: preliminary results from a multicenter randomized controlled clinical trial**
Journal of Clinical Periodontology, 2017 Feb;44(2):204-214
- **E)** Crespi R et al. - **Corticocancellous porcine bone in the healing of human extraction sockets: combining histomorphometry with osteoblast gene expression profiles in vivo**
Int Journal of Oral and Maxillofacial Implants, 2011 Jul - Aug; 26(4):866-72
- **F)** Crespi R et al. - **Comparison of magnesium-enriched hydroxyapatite and porcine bone in human extraction socket healing: a histologic and histomorphometric evaluation**
Int Journal of Oral and Maxillofacial Implants, 2011 Sep-Oct;26(5):1057-62
- **G)** Barone A et al. - **Xenograft versus extraction alone for ridge preservation after tooth removal: a clinical and histomorphometric study**
Journal of Periodontology, 2008 Aug; 79(8):1370-7

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- **A)** Cassetta M et al. - **Bone formation in sinus augmentation procedures using autologous bone, porcine bone, and a 50 : 50 mixture: a human clinical and histological evaluation at 2 months**
Clinical Oral Implants Research, 2015 Oct;26(10):1180-4
- **B)** Scarano A et al. - **Maxillary sinus augmentation in humans using cortical porcine bone: a histological and histomorphometrical evaluation after 4 and 6 months**
Clinical Implant Dentistry and Related Research, 2011 Mar; 13(1):13-18
- **C)** Orsini G et al. - **Histologic and ultrastructural analysis of regenerated bone in maxillary sinus augmentation using a porcine bone-derived biomaterial**
Journal of Periodontology, 2006 Dec; 77(12):1984-90
- **D)** Barone A et al. - **A 6-month histological analysis on maxillary sinus augmentation with and without use of collagen membranes over the osteotomy window: randomized clinical trial**
Clinical Oral Implants Research, 2013 Jan; 24(1):1-6
- **E)** Iezzi G et al. - **Comparative histological results of different biomaterials used in sinus augmentation procedures: a human study at 6 months**
Clinical Oral Implants Research, 2012 Dec;23(12):1369-76
- **F)** Silvestri M et al. - **Simultaneous sinus augmentation with implant placement: histomorphometric comparison of two different grafting materials. A multicenter double-blind prospective randomized controlled clinical trial**
Int Journal of Oral and Maxillofacial Implants, 2013 Mar-Apr; 28(2):543-9
- **G)** Barone A et al. - **Maxillary sinus augmentation using prehydrated corticocancellous porcine bone: histomorphometric evaluation after 6 months**
Clinical Implant Dentistry and Related Research, 2012 Jun;14(3):373-9

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
BIBLIOGRAFIA

- (1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (2) KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
BMC ORAL HEALTH, 2017 OCT 5;17(1):126
- (3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (5) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVEST, 2017 2018 JAN;22(1):505-513
- (6) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6
- (7) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (8) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- (9) BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY
CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112
- (10) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

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



(3) TALLARICO M, KHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM
SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):263-275


(4) FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESSENCE INT, 2014 NOV-DEC;45(10):853-60

(5) FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES
QUINTESSENCE INT, 2013;44(3):243-6

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
BIBLIOGRAFIA

- (1) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (2) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58
- (3) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT AND RELAT RES, 2013 OCT;15(5):707-13
- (4) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (5) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85
- (6) HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197
- (7) ROSSI R, FOCE E, SCOLAVINO S
THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT
J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41

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INNOVAZIONE

A close-up photograph of a scientist in a white lab coat, hairnet, and face mask, looking through a microscope. The scientist is positioned on the right side of the frame, with their face partially obscured by the mask and hairnet. The microscope is the central focus, with its eyepieces and objective lenses clearly visible. The background is a blurred laboratory setting with various pieces of equipment. The word "INNOVAZIONE" is overlaid in large, bold, blue capital letters on the left side of the image.

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BIBLIOGRAFIA

- (1) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (2) BARONE A, ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS
CURR PHARM BIOTECHNOL, 2017;18(1):64-75
- (3) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (4) BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73
- (7) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (8) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (9) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (10) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4

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BIBLIOGRAFIA

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

(3) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(4) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6

(5) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

(6) BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9

(7) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

(8) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

(9) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42

(10) MJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR MATERIALS, 2017 SEP 8;10(9)

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BIBLIOGRAFIA

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H

PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE

J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A

HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL

J PERIODONTOL, 2006 DEC;77(12):1984-90

(3) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C

ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION

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VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT

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MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS

CURR PHARM BIOTECHNOL, 2017,18,33-44

(6) MIZUNO M, FUJISAWA R, KUBOKI Y

TYPE I COLLAGEN-INDUCED OSTEOBLASTIC DIFFERENTIATION OF BONE-MARROW CELLS MEDIATED BY COLLAGEN-A2B1 INTEGRIN INTERACTION

J CELL PHYSIOL, 2000 AUG;184(2):207-13

(7) HSU FY, CHUEH SC, WANG YJ

MICROSPHERES OF HYDROXYAPATITE/RECONSTITUTED COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH

BIOMATERIALS 1999, 20:1931-1936

(8) ABDELGAWAD ME, SØE K, ANDERSEN TL, MERRILL DM, CHRISTIANSEN P, KJÆRGAARD-ANDERSEN P, DELAISSE JM

DOES COLLAGEN TRIGGER THE RECRUITMENT OF OSTEOBLASTS INTO VACATED BONE RESORPTION LACUNAE DURING BONE REMODELING?

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BIBLIOGRAFIA

(1) GRIFFITH LG, NAUGHTON G
TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES
SCIENCE 2002, 295:1009-14

(2) REDDI AH
MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS
TISSUE ENG 2000, 6(4):351-59

(3) NAKASHIMA N, REDDI AH
THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING
NAT BIOTECHNOL 2003, 9:1025-32

(4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE
ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS
J BIOMED BIOTECHNOL 2004, 1:24-34

(5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL[®] INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

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(2) NANMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

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CASSETTA M, CALASSO S, VOZZA I, DELLAQUILA D
**REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH
CYLINDRICAL SANDBLASTED AND ACID ETCHED
IMPLANTS: A PILOT STUDY**
EUR J IMPLANT PROSTHODONTICS, 2005;3(1):133-144

ARCURI C, CECCHETTI F, GERMANO F, MOTTA A,
SANTACROCE C
**CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC
BONE SUBSTITUTE USED AS A FILLER IN
POSTEXTRACTIVE ALVEOLUS**
MINERVA STOMATOL, 2005 JUN;54(6):351-62

BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R,
COVANI U
**MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND
HISTOMORPHOMETRIC ANALYSIS**
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG;
20(4):519-25

RINNA C, UNGARI C, SALTAREL A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
J CRANIOFAC SURG, 2005 NOV; 16(6):968-72

BARONE A, AMERI S, COVANI U
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT
OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE
ANALYSIS**
EUR J IMPLANT PROSTHODONTICS, 2006;2: 99-106

BARONE A, SANTINI S, SBORDONE L, CRESPI R, COVANI U
**A CLINICAL STUDY OF THE OUTCOMES AND
COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS
AUGMENTATION**
INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5

COVANI U, BARONE A, CORNELINI R, CRESPI R
**CLINICAL OUTCOME OF IMPLANTS PLACED
IMMEDIATELY AFTER IMPLANT REMOVAL**
J PERIODONTOL, 2006 APR;77(4):722-7

ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S,
PIATTELLI A
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF
REGENERATED BONE IN MAXILLARY SINUS
AUGMENTATION USING A PORCINE BONE-DERIVED
BIOMATERIAL**
J PERIODONTOL, 2006 DEC;77(12):1984-90

TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D,
D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
**THE PERFORMANCE OF HUMAN PERIODONTAL
LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC
BIOMATERIALS**
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20
(1 SUPPL 1):87-91

BARONE A, COVANI U
**MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH
NON- VASCULARIZED AUTOGENOUS BLOCK BONE:
CLINICAL RESULTS**
J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46

DEL CORSO M
**SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN:
CLINICAL EVIDENCES**
COSMETIC DENTISTRY, 2008, 3:16-20

BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L,
GHERLONE E, COVANI U
**OSTEOTOMY AND MEMBRANE ELEVATION DURING THE
MAXILLARY SINUS AUGMENTATION PROCEDURE. A
COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS.
CONVENTIONAL ROTATIVE INSTRUMENTS**
CLIN ORAL IMPLANTS RES, 2008 MAY;19(5):511-5

BARONE A, CORNELINI R, CIAGLIA R, COVANI U
**IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS
AND SIMULTANEOUS OSTEOTOME SINUS FLOOR
ELEVATION: A CASE SERIES**
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN;
28(3):283-9

BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO
GUIRADO JL, COVANI U
**XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE
PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL
AND HISTOMORPHOMETRIC STUDY**
J PERIODONTOL, 2008 AUG;79(8):1370-7

COVANI U, CORNELINI R, BARONE A
**BUCCAL BONE AUGMENTATION AROUND IMMEDIATE
IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A
MODIFIED APPROACH**
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT;
23(5):841-6

CARDAROPOLI D, CARDAROPOLI G
**PRESERVATION OF THE POSTEXTRACTION ALVEOLAR
RIDGE: A CLINICAL AND HISTOLOGIC STUDY**
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT;
28(5):469-77

NANNMARK U, SENNERBY L
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND
COLLAGENATED CORTICO-CANCELLOUS PORCINE
BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V
**REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED
INTO THE MAXILLARY SINUS. A CLINICAL AND
HISTOLOGIC CASE REPORT**
J OSSEOINTEGR, 2009;1(1):35-40

COVANI U, MARCONCINI S, CRESPI R, BARONE A
**IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A
FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE
REPORT**
J ORAL IMPLANTOL, 2009; 35(4):189-95

FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F,
FIGUEIREDO H
**PHYSICO-CHEMICAL CHARACTERIZATION OF
BIOMATERIALS COMMONLY USED IN DENTISTRY AS
BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**
J BIOMED MATER RES B APPL BIOMATER, 2010FEB;
92(2):409-19

GRENGA PL, REALE G, COFONE C, MEDURI A, CERUTI P,
GRENGA R
**HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30
PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL
BLOW-OUT FRACTURE**
OPHTHAL PLAST RECONSTR SURG, 2009 MAR-APR;
25(2):123-5

CRESPI R, CAPPARÈ P, GHERLONE E
**DENTAL IMPLANTS PLACED IN EXTRACTION SITES
GRAFTED WITH DIFFERENT BONE SUBSTITUTES:
RADIOGRAPHIC EVALUATION AT 24 MONTHS**
J PERIODONTOL, 2009 OCT; 80(10):1616-1621

RINNA C, REALE G, FORESTA E, MUSTAZZA MC
**MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE
BONE CORTEX**
J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4

CARDAROPOLI D, CARDAROPOLI G
**HEALING OF GINGIVAL RECESSIONS USING A
COLLAGEN MEMBRANE WITH A THE MINERALIZED
XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL
TRIAL**
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB;
29(1):59-67

NANNMARK U, AZARMEHR I
**SHORT COMMUNICATION: COLLAGENATED CORTICO-
CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN
RABBIT MAXILLARY DEFECTS**

CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A,
PERROTTI V, PIATTELLI M, IEZZI G
**PORCINE BONE USED IN SINUS AUGMENTATION
PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL
EVALUATION**
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N
**PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A
NEW IMPLANT 3D NAVIGATION SYSTEM**
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88

SCARANO A, CARINCI F, ASSENZA B, PIATTELLI M, MURMURA
G, PIATTELLI A
**VERTICAL RIDGE AUGMENTATION OF ATROPHIC
POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE
WITH A XENOGRAFT WITHOUT MINISCREWS AND
MINIPLATES: CASE SERIES**

CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30

PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A,
VERROCCHI D, VOLPE S, SENNERBY L
**A COLLAGENATED PORCINE BONE SUBSTITUTE FOR
AUGMENTATION AT NEOS IMPLANT SITES: A
PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY
WITH HISTOLOGY**

CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58

SANTAGATA M, GUARINIELLO L, TARTARO G
**A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE)
TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS.
A CASE REPORT**
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N°:114-9

BARONE A, RICCI M, CALVO GUIRADO JL, COVANI U
**BONE REMODELLING AFTER REGENERATIVE
PROCEDURES AROUND IMPLANTS PLACED IN FRESH
EXTRACTION SOCKETS: AN EXPERIMENTAL STUDY IN
BEAGLE DOGS**
CLIN ORAL IMPLANTS RES, 2011 OCT; 22(10):1131-7

SCARANO A, PIATTELLI A, PERROTTI V, MANZON I, IEZZI G
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING
CORTICAL PORCINE BONE: A HISTOLOGICAL AND
HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6
MONTHS**
CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18

CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E,
BENASCIOTTI E, GHERLONE E
**CORTICOCANCELLOUS PORCINE BONE IN THE
HEALING OF HUMAN EXTRACTION SOCKETS:
COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST
GENE EXPRESSION PROFILES IN VIVO**
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG;
26(4):866-72

HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
**ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION
WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS
ELEVATION TECHNIQUE". A CASE-COHORT STUDY**
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A,
SHIBLI JA, PERROTTI V
**COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT
BIOMATERIALS USED IN SINUS AUGMENTATION
PROCEDURES: A HUMAN STUDY AT 6 MONTHS**
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76

SLOTTE C, LINDFORS N, NANNMARK U
**SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE
DEFECTS WITH PREHYDRATED AND COLLAGENATED
PORCINE BONE AND COLLAGEN BARRIERS: CASE
PRESENTATIONS**
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23

BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA
A, COVANI U
**A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY
SINUS AUGMENTATION WITH AND WITHOUT USE OF
COLLAGEN MEMBRANES OVER THE OSTEOTOMY
WINDOW: RANDOMIZED CLINICAL TRIAL**
CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6

SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
**IMMEDIATE LOADING OF DENTAL IMPLANT AFTER
SINUS FLOOR ELEVATION WITH OSTEOTOME
TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY
RADIOGRAPHIC RESULTS**
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A
SOFT CORTICAL MEMBRANE VERSUS EXTRACTION
ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL
STUDY IN HUMANS**
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13

RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ
SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA
DORADO C
**ULTRASTRUCTURAL STUDY BY BACKSCATTERED
ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS
OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL
DEGRADATION OF PORCINE XENOGRAFTS USED IN
MAXILLARY SINUS FLOOR ELEVATION**
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

CASSETTA M, RICCI I, IEZZI G, DELLAQUILA D, PIATTELLI A,
PERROTTI V
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS
INSERTED WITH A SIMULTANEOUS GRAFTING
PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**
INT J PERIODONTICS RESTORATIVE DENT, 2012
OCT;32(5):581-9

BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U
A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION
J PERIODONTOL, 2012 JUL;83(7):836-46

ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33

FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
IMPLANT DENT, 2012 JUN;21(3):190-5

BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8

BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (S), 103-107

FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES
QUINTESSENCE INT, 2013;44(3):243-6

SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9

WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

RODRIGUEZ JG, ELDIBANY RM
VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION
INT J ORAL MAXILLOFAC SURG, 2013 SEP;42(9):1060-6
FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1;33(6):3506-13

FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61

TRAINI T, PIATTELLI A, CAPUTI S, DEGIDI M, MANGANO C, SCARANO A, PERROTTI V, IEZZI G
REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS
CLIN IMPLANT DENT RELAT RES, 2015 FEB;17(1):150-62

KOLMAS J, SZWAJA M, KOLODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENIC BIOMATERIALS USED AS BONE SUBSTITUTES
J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41

PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARAUSSÉ C, BUTI J, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2013 WINTER;6(4):343-357

FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KESCHULL M
PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS - AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):799-805

BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

BARONE A, RICCI M, ROMANOS GE, TONELLI P, ALFONSI F, COVANI U
BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):823-30

BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26

CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4

GHEÑO E, PALERMO A, BUFFOLI B, RO DELLA LF
THE EFFECTIVENESS OF THE USE OF XENOGENIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES
J OSSEOINTEGRATION 2014;6(2):37-42

FALISI G, GALLI M, VITTORINI-VELASQUEZ P, GALLEGOS-RIVERA JC, MINASI R, DE BIASE A, DI PAOLO C
USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE
ACTA ODONTOL LATINOAM 2013;26(3):167-172

FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: TWO CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESSENCE INT, 2014 NOV-DEC; 45(10):853-60

FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42

SCARANO A, MURMURA G, SINJARI B, ASSENZA B, SOLLAZZO V, SPINELLI G, CARINCI F
EXPANSION OF THE ALVEOLAR BONE CREST WITH ULTRASONIC SURGERY DEVICE: CLINICAL STUDY IN MANDIBLE
INT J IMMUNOPATHOL PHARMACOL, 2011 APR-JUN; 24(2 SUPPL):71-5

SCARANO A, PIATTELLI A, MURMURA G, IEZZI G, ASSENZA B, MANCINO C
DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGIC CASE SERIES
INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9

LORENZON G, BUTTARELLO GM, CHESSA G
CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION
DENTISTRY, 2015, 5:2

THALMAIR T, FICKL S, SCHNEIDER D, HINZE M, WACHTEL H
DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES - A VOLUMETRIC STUDY
J CLIN PERIODONTOL, 2013 JUL;40(7):721-7

MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATTELLI A, TRUBIANI O
OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON
J PERIODONTAL RES, 2016 FEB;51(1):112-24

SCARANO A, PIATTELLI A, IEZZI G, VARVARA G
SPONTANEOUS BONE FORMATION ON THE MAXILLARY SINUS FLOOR IN ASSOCIATION WITH SURGERY TO REMOVE A MIGRATED DENTAL IMPLANT: A CASE REPORT
MINERVA STOMATOL, 2014 OCT;63(10):351-9

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2016 NOV;27(11):E105-E115

ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244

OZEL B, FINDIKCIOGLU K, SEZGIN B, GUNAY K, BARUT I, OZMEN S
A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE
J CRANIOMAXILLOFAC SURG, 2015 OCT;43(8):1583-8

CORBELLA S, TASCHEI S, WEINSTEIN R, DEL FABBRO M
HISTOMORPHOMETRIC OUTCOMES AFTER LATERAL SINUS FLOOR ELEVATION PROCEDURE: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS
CLIN ORAL IMPLANTS RES, 2016 SEP;27(9):1106-22

BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY
CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112

CAKIR M, KARACA IR, AYŞEGÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION
CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F
REGENERATION OF ATROPHIC CRISTAL RIDGES WITH RESORBABLE LAMINA: TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS 2015 JUL-SEP;29(3 SUPPL 1):97-100

ETTORRE V, DE MARCO P, ZARA S, PERROTTI V, SCARANO A, DI CRESCENZO A, PETRINI M, HADAD C, BOSCO D, ZAVAN B, VALBONETTI L, SPOTO G, IEZZI G, PIATTELLI A, CATALDI A, FONTANA A
IN VITRO AND IN VIVO CHARACTERIZATION OF GRAPHENE OXIDE COATED PORCINE BONE GRANULES CARBON, JULY 2016, VOLUME 103, PAGES 291-298

ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68

SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870

FELICE P, ZUCCHELLI G, CANNIZZARO G, BARAUSSE C, DIAZZI M, TRULLENQUE-ERIKSSON A, ESPOSITO M
IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 4-MONTH POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):233-247

ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

BARONE A, TOTI P, MARCONCINI S, DERCHI G, MARCHIONNI S, COVANI U
ESTHETIC OUTCOME OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS BY CLINICIANS WITH OR WITHOUT EXPERIENCE: A MEDIUM-TERM RETROSPECTIVE EVALUATION
INT J ORAL MAXILLOFAC IMPLANTS, 2016;31(6)

TALLARICO M, KHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM
SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):263-275

LOPEZ MA, MANZULLI N, CASALE M, ORMIANER Z, CARINCI F
THE USE OF RESORBABLE HETEROLOGOUS CORTICAL LAMINA AS A NEW SINUS LIFT FLOOR: A TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):75-79

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

ESPOSITO M, ZUCCHELLI G, BARAUSSE C, PISTILLI R, TRULLENQUE-ERIKSSON A, FELICE P
FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 4-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(4):393-409

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U
CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAPTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL
J CLIN PERIODONTOL, 2017 FEB;44(2):204-214

BARONE A, ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS
CURR PHARM BIOTECHNOL, 2017;18(1):64-75

BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

BARONE A, TOTI P, FUNEL N, CAMPANI D, COVANI U
EXPRESSION OF SP7, RUNX1, DLX5, AND CTNBN1 IN HUMAN MESENCHYMAL STEM CELLS CULTURED ON XENOGENIC BONE SUBSTITUTE AS COMPARED WITH MACHINED TITANIUM
IMPLANT DENT, 2014 AUG;23(4):407-15

ESPOSITO M, ZUCCHELLI G, CANNIZZARO G, CHECCHI L, BARAUSSE C, TRULLENQUE-ERIKSSON, FELICE P
IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 1-YEAR POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(1):11-26

SCARANO A, CRINCOLI V, DI BENEDETTO A, COZZOLINO V, LORUSSO F, PODALIRI VULPIANI M, GRANO M, KALEMAJ Z, MORI G, GRASSI FR
BONE REGENERATION INDUCED BY BONE PORCINE BLOCK WITH BONE MARROW STROMAL STEM CELLS IN A MINIPIG MODEL OF MANDIBULAR "CRITICAL SIZE" DEFECT STEM CELLS INT, 2017;2017:9082869

SCARANO A
TRADITIONAL POSTEXTRACTIVE IMPLANT SITE PREPARATION COMPARED WITH PRE-EXTRACTIVE INTERRADICULAR IMPLANT BED PREPARATION IN THE MANDIBULAR MOLAR REGION, USING AN ULTRASONIC DEVICE: A RANDOMIZED PILOT STUDY
INT J ORAL MAXILLOFAC IMPLANTS, 2017 MAY/JUN;32(3):655-660

BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLINICAL ORAL INVESTIGATIONS, 2018 JAN;22(1):505-513

IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44

FELICE P, BARAUSSE C, BARONE A, ZUCCHELLI G, PIATTELLI M, PISTILLI R, IPPOLITO DR, SIMION M
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

IIDA T, CARNEIRO MARTINS NETO E, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP
INFLUENCE OF A COLLAGEN MEMBRANE POSITIONED SUBJACENT THE SINUS MUCOSA FOLLOWING THE ELEVATION OF THE MAXILLARY SINUS. A HISTOMORPHOMETRIC STUDY IN RABBITS
CLIN IMPLANT DENT RELAT RES, 2017 JUN 7, EPUB AHEAD OF PRINT

DE MARCO P, ZARA S, DE COLLI M, RADUNOVIC M, LAZOVIC V, ETORRE V, DI CRESCENZO A, PIATTELLI A, CATALDI A, FONTANA A
GRAPHENE OXIDE IMPROVES THE BIOCOMPATIBILITY OF COLLAGEN MEMBRANES IN AN IN VITRO MODEL OF HUMAN PRIMARY GINGIVAL FIBROBLASTS
BIOMED MATER, 2017 SEP 13;12(5):055005

MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR
MATERIALS, 2017 SEP 8;10(9)

ROSSI R, FOCE E, SCOLAVINO S
THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT
J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41

CHECCHI V, FELICE P, ZUCCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278

CRESPI R, CAPPARÈ P, GHERLONE E
COMPARISON OF MAGNESIUM-ENRICHED HYDROXYAPATITE AND PORCINE BONE IN HUMAN EXTRACTION SOCKET HEALING: A HISTOLOGIC AND HISTOMORPHOMETRIC EVALUATION
INT J ORAL MAXILLOFAC IMPLANTS, 2011 SEP-OCT;26(5):1057-62

CORBELLA S, TASCHIERI S, FRANCIETTI L, WEINSTEIN R, DEL FABBRO M
HISTOMORPHOMETRIC RESULTS AFTER POSTEXTRACTION SOCKET HEALING WITH DIFFERENT BIOMATERIALS: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2017 SEP/OCT;32(5):1001-1017

RADUNOVIC M, DE COLLI M, DE MARCO P, DI NISIO C, FONTANA A, PIATTELLI A, CATALDI A, ZARA S
GRAPHENE OXIDE ENRICHMENT OF COLLAGEN MEMBRANES IMPROVES DPSCS DIFFERENTIATION AND CONTROLS INFLAMMATION OCCURRENCE
J BIOMED MATER RES A, 2017 AUG;105(8):2312-2320

KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
BMC ORAL HEALTH, 2017 OCT 5;17(1):126

TROIANO G, ZHURAKIVSKA K, LO MUZIO L, LAINO L, CICCÌ M, LO RUSSO L
COMBINATION OF BONE GRAFT AND RESORBABLE MEMBRANE FOR ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW, META-ANALYSIS AND TRIAL SEQUENTIAL ANALYSIS
J PERIODONTOL, 2017 SEP 12:1-17. EPUB AHEAD OF PRINT

ROSSI R, LONGO E, MIJIRITSKY E
A NEW INTERPRETATION OF GUIDED IMPLANT SURGERY TO ACHIEVE AN OPTIMAL RESULT IN THE ESTHETIC ZONES
MEDICAL RESEARCH ARCHIVES, 2017 APRIL, VOL. 5, ISSUE 4

BOLLE C, FELICE P, BARAUSSE C, PISTILLI V, TRULLENQUE-ERIKSSON A, ESPOSITO M
FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):31-47

ESPOSITO M, DAVÒ R, MARTI PAGES C, FERRER FUENTES A, BARAUSSE C, PISTILLI R, IPPOLITO DR, FELICE P
IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 4 MONTHS POST-LOADING RESULTS FROM A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):11-28

GASTALDI G, FELICE P, PISTILLI V, BARAUSSE C, IPPOLITO DR, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 x 5 MM IMPLANTS WITH A NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):49-61

DIOMEDE F, D'AUROA M, GUGLIANDOLO A, MERCIARO I, ORSINI T, GATTA V, PIATTELLI A, TRUBIANI O, MAZZON E
BIOFUNCTIONALIZED SCAFFOLD IN BONE TISSUE REPAIR
INT J OF MOLECULAR SCIENCES, 2018, 19, 1022

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


Tecross[®] s.r.l.

Via Nurivalle, 8
10094 Giaveno (TO), Italy
Tel./Fax. +39 011 976 6684

info@tecross.com
www.tecross.com
BIOMATERIALS ENGINEERING

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All rights reserved.